

MULTIMODAL PNEUMONIA DETECTION SYSTEM

Complete Technical Report — Dual-Branch Fusion Pipeline

CSV Clinical Data · Chest X-Ray Images · Fusion Architecture · Full ML Pipeline

University of Saida, Algeria · AI & Medical Imaging · 2025–2026

Abstract. This report describes the complete design and implementation plan for a multimodal AI system that combines two independent data sources — a structured clinical CSV dataset (patient vitals and symptoms) and a chest X-ray image dataset — to detect and classify pediatric pneumonia. The report covers: (1) both dataset descriptions and download links, (2) the pediatric filtering problem in the CSV data and its solution, (3) the full dual-branch parallel training pipeline, (4) input/output specifications at every stage, (5) all fusion strategies with code, (6) the complete tabular ML pipeline with 7 models and cross-validation, and (7) implementation roadmap. Current project status: literature review and preprocessing complete — model training is the next step.

Branch	Dataset Name	Type	Kaggle Link
1 — Clinical	essienmary/pneumonia-dataset	CSV — vitals, symptoms, age, label	kaggle.com/datasets/essienmary/pneumonia-dataset
2 — X-Ray	paultimothymooney/chest-xray-pneumonia	JPEG chest X-rays — Normal / Pneumonia	kaggle.com/datasets/paultimothymooney/chest-xray-pneumonia

PROJECT TEAM — University of Saida, Algeria

Role	Name	Contact / Institution
Data Engineer	Bouhmidi Amina Maroua	University of Saida
ML Engineer	Labani Nabila Nour Elhouda	University of Saida
Project Manager	Kassouar Fatima	kassouarfatima48@gmail.com
Business Model	Miloudi Meroua Amira	University of Saida
Academic Supervisor	Dr. Abderrahmane Khat	University of Saida, Algeria
Medical Advisor	Dr. Aimer Mohammed Djamel Eddine	

TABLE OF CONTENTS	
Section	Title
1	Dataset Descriptions — CSV Clinical + X-Ray Images
2	The Problem: CSV Data Is Not Pediatric-Only + Solution
3	Full Multimodal Pipeline — Big Picture Overview
4	Preprocessing Pipeline — Branch A (CSV) + Branch B (X-Ray)
5	Branch A — Tabular ML: 7 Models + Cross-Validation + Code
6	Branch B — CNN Image Models: 8 Architectures + Transfer Learning
7	Fusion Strategies — Late Fusion + Intermediate Fusion + Code
8	Evaluation Metrics + Validation Strategy
9	Implementation Roadmap + Current Project Status
10	Conclusion + Recommendations

SECTION 1 — DATASET DESCRIPTIONS

1.1 Branch A — Clinical CSV Dataset

Name: Pneumonia Patient Dataset | **Author:** essienmary

Link: <https://www.kaggle.com/datasets/essienmary/pneumonia-dataset>

Type: Structured tabular CSV — patient vitals, clinical symptoms, lab indicators, diagnosis label

This dataset contains structured clinical records for patients with and without pneumonia. Each row represents one patient with measurements of vital signs and symptoms collected at the time of hospital visit. The dataset includes patients of **all age groups** — from infants to elderly — which creates a critical problem for this project (see Section 2).

Column Name	Data Type	Description	Clinical Role
patient_id	Integer / ID	Unique patient identifier	Drop — not a predictive feature
age	Numeric (years)	Patient age	Critical for pediatric filter
sex	Categorical M/F	Biological sex	Encode: M=0, F=1
temperature	Numeric (°C)	Body temperature	Fever >38.5°C = active infection
spo2	Numeric (%)	Blood oxygen saturation	<95% = respiratory compromise
heart_rate	Numeric (bpm)	Heart beats per minute	Tachycardia common in pneumonia
respiratory_rate	Numeric (/min)	Breathing rate	>40/min in children = tachypnea
cough	Binary 0/1	Cough present or absent	Primary pneumonia symptom
label	Binary 0/1	0=Normal, 1=Pneumonia	Main classification target
pneumonia_severity	Ordinal 4-class	none / mild / moderate / severe	Secondary target for severity scoring

1.2 Branch B — Chest X-Ray Image Dataset

Name: Chest X-Ray Images (Pneumonia) — Kermay et al., 2018

Link: <https://www.kaggle.com/paultimothymooney/chest-xray-pneumonia>

Type: JPEG chest X-ray images — 5,856 images — **Pediatric patients aged 1–5 years only**

This dataset contains frontal chest X-rays collected from the Guangzhou Women and Children's Medical Center, China. All images were graded by two expert physicians and a third for verification. Images are stored in train/val/test folders with NORMAL and PNEUMONIA subfolders.

Split	NORMAL	PNEUMONIA	Total	Problem / Note
Train	1,341	3,875	5,216	Class imbalance 1:2.89 — needs class weighting
Validation	8	8	16	TOO SMALL — must be replaced with proper stratified split
Test	234	390	624	Final held-out evaluation set — do not touch during training
TOTAL	1,583	4,273	5,856	Stratified 70/15/15 split recommended (team fix)

SECTION 2 — THE PROBLEM: CSV IS NOT PEDIATRIC-ONLY + SOLUTION

2.1 The Core Problem

CRITICAL MISMATCH: The X-ray dataset (Branch B) contains ONLY pediatric patients aged 1–5 years. The CSV dataset (Branch A) contains patients of ALL ages — adults, elderly, and children mixed together. Using the full CSV with pediatric X-rays makes the fusion model clinically inconsistent and scientifically invalid.

Pneumonia has completely different clinical presentations depending on age. Vital sign ranges, typical pathogens, and X-ray patterns differ fundamentally between a 2-year-old child and a 70-year-old adult. A model trained on mixed-age clinical data and pediatric X-rays would learn incorrect cross-modal relationships.

Age Group	Normal SpO2	Resp. Rate (/min)	Typical Pathogen	X-Ray Pattern
Infant 0–1 yr	95–100%	30–60	RSV, viral	Bilateral haziness
Child 1–5 yr ★	97–100%	24–40	Viral, S. pneumoniae	Lobar / perihilar
Adult 18–60 yr	95–99%	12–20	S. pneumoniae, atypical	Lobar consolidation
Elderly 60+ yr	92–97%	15–25	Multi-resistant strains	Diffuse bilateral

★ Only this age group exists in the X-ray dataset — the two branches must match this population.

2.2 The Solution — Filter CSV to Pediatric Ages 1–5

SOLUTION: Before any training, filter the CSV dataset to keep ONLY rows where age is between 1 and 5 years. This aligns both branches on the same patient population and makes the fusion clinically valid. If too few samples remain after filtering, expand to ages 0–12.

Step-by-Step Filtering Code:

```
import pandas as pd # Step 1 - Load the full CSV df = pd.read_csv('pneumonia_dataset.csv') print(f'Full dataset shape: {df.shape}') print(f'Age range: {df["age"].min():.0f} to {df["age"].max():.0f} years') print(f'Age distribution:\n{df["age"].describe()}') # Step 2 - Check class balance BEFORE filter print(f'\nLabel balance (full):') print(df['label'].value_counts()) # Step 3 - Apply pediatric filter (ages 1-5, same as X-ray dataset) df_kids = df[(df['age'] >= 1) & (df['age'] <= 5)].copy() print(f'\nPediatric dataset (ages 1-5): {len(df_kids)} rows') print(f'Label balance (pediatric):') print(df_kids['label'].value_counts()) # Step 4 - If sample size is too small, expand to all pediatric (0-12) if len(df_kids) < 200: print('WARNING: too few samples. Expanding to ages 0-12.') df_kids = df[(df['age'] >= 0) & (df['age'] <= 12)].copy() print(f'Expanded pediatric dataset: {len(df_kids)} rows') # Step 5 - Save filtered version for use in Branch A df_kids.to_csv('pneumonia_dataset_pediatric.csv', index=False) print('Saved: pneumonia_dataset_pediatric.csv')
```

Step	Action	Why
1	Load full CSV and check age distribution	Understand the age spread in the dataset
2	Filter: age >= 1 AND age <= 5	Match the X-ray dataset population exactly
3	Check class balance after filter	Ensure enough positive and negative cases remain
4	Expand to 0–12 if fewer than 200 rows remain	Statistical minimum for reliable model training
5	Save as pneumonia_dataset_pediatric.csv	Use this file exclusively for Branch A training

SECTION 3 — FULL MULTIMODAL PIPELINE OVERVIEW

3.1 The Big Picture — Mimicking a Real Doctor

A real physician diagnoses pneumonia by consulting two sources simultaneously: the **patient chart** (vitals, symptoms, history) AND the **chest X-ray**. This system mirrors that workflow exactly. Two AI models are trained in parallel, each on one data source, and their outputs are fused for a final prediction. This is called a **Multimodal Fusion Architecture**.




Figure 1: Complete multimodal pipeline. Both branches trained in parallel on independent datasets, fused for final prediction.

3.2 Input/Output Specification at Each Stage

Stage	Branch	Input	Output	Shape
1. Raw CSV	A	pneumonia_dataset.csv	Full DataFrame — all ages	N rows × 10 cols
2. Filter	A	All ages mixed	Pediatric only (ages 1–5)	n rows × 10 cols
3. Preprocess A	A	Raw pediatric rows	Normalized + encoded feature matrix	[n, 16] float32
4. Raw X-Ray	B	JPEG files in folders	Pixel arrays	[H, W, 3] uint8
5. Preprocess B	B	Raw JPEG	Normalized float image	[224, 224, 3] float32
6. Train Tabular	A	Feature matrix [n, 16]	7 trained models + CV metrics	prob p1 ∈ [0,1]
7. Train CNN	B	Image batch [B, 224, 224, 3]	8 trained CNNs + metrics	prob p2 ∈ [0,1]
8. Select Best	A+B	CV and val AUC scores	1 tabular model + 1 CNN selected	Two models
9. Late Fusion	A+B	p1 (tabular) + p2 (CNN)	Weighted average p_final	scalar ∈ [0,1]
10. Decision	A+B	p_final	PNEUMONIA/NORMAL + severity	Label + score

SECTION 4 — PREPROCESSING PIPELINES

4.1 Branch A — Clinical CSV Preprocessing

All preprocessing steps are applied **inside a scikit-learn Pipeline** to prevent data leakage. SMOTE is applied only inside each cross-validation fold on the training portion. The scaler is fit on training data only and applied to validation/test.

Step	Operation	Detail / Reason
1	Load pediatric CSV	Use pneumonia_dataset_pediatric.csv (ages 1–5 only)
2	Drop patient_id	Not a predictive feature — unique identifier only
3	Encode sex	LabelEncoder: Male=0, Female=1
4	Encode severity	OrdinalEncoder: none=0, mild=1, moderate=2, severe=3
5	Median imputation	SimpleImputer(strategy=median) for missing vitals
6	Feature engineering	Create: danger_score = (fever) + (low_spo2) + (tachypnea)
7	Stratified 80/20 split	Preserve class ratio in train and test
8	StandardScaler	Fit on X_train only, transform X_train and X_test
9	SMOTE (train fold only)	Synthetic oversampling — applied inside CV folds only to prevent leakage

Full Preprocessing Code — Branch A:

```
import pandas as pd
import numpy as np
from sklearn.preprocessing import LabelEncoder, OrdinalEncoder, StandardScaler
from sklearn.impute import SimpleImputer
from sklearn.model_selection import train_test_split
from imblearn.over_sampling import SMOTE
from imblearn.pipeline import Pipeline as ImbPipeline

# Load filtered pediatric CSV
df = pd.read_csv('pneumonia_dataset_pediatric.csv')

# Encode categorical columns
le = LabelEncoder()
df['sex'] = le.fit_transform(df['sex'])
oe = OrdinalEncoder(categories=[['none', 'mild', 'moderate', 'severe']])
df['pneumonia_severity'] = oe.fit_transform(df[['pneumonia_severity']])

# Feature engineering
df['high_fever'] = (df['temperature'] > 38.5).astype(int)
df['low_spo2'] = (df['spo2'] < 95).astype(int)
df['tachypnea'] = (df['respiratory_rate'] > 40).astype(int)
df['danger_score'] = df['high_fever'] + df['low_spo2'] + df['tachypnea']

# Drop ID and prepare X, y
df.drop('patient_id', axis=1, inplace=True, errors='ignore')
X = df.drop('label', axis=1)
y = df['label']

# Stratified 80/20 split
X_train, X_test, y_train, y_test = train_test_split(X, y, test_size=0.2, random_state=42, stratify=y)

print(f'Train: {len(X_train)} | Test: {len(X_test)}')
print(f'Features: {X_train.shape[1]}')
print(f'Class balance (train): {y_train.value_counts().to_dict()}')
```

4.2 Branch B — X-Ray Image Preprocessing

Step	Operation	Value / Detail
1	Fix validation split	Replace 16-image val set with stratified 70/15/15 (4099/878/879)
2	Resize	224×224 px — standard ImageNet input size
3	Convert to float32	Required — float16 causes sigmoid overflow (NaN loss)
4	Normalize	÷255 → [0,1] for all models except MobileNetV2 (÷127.5-1 → [-1,1])
5	Augmentation (train only)	Horizontal flip, rotate ±36°, zoom ±15%, brightness ±15%
6	Class weighting	w_normal=2.89, w_pneumonia=1.0 to compensate 1:2.89 imbalance
7	Batching	Batch size 32 (GPU) or 16 (CPU). Prefetch for speed.

CRITICAL: MobileNetV2 requires pixel values in $[-1, 1]$ (divide by 127.5 then subtract 1). All other models use $[0, 1]$ (divide by 255). Mixing these reduces accuracy by 5–15% silently. Also: the final Dense and Sigmoid layers **MUST** use `dtype=float32` to avoid NaN loss under FP16.

SECTION 5 — BRANCH A: 7 TABULAR MODELS WITH CROSS-VALIDATION

5.1 Model Overview

Each of the 7 models is wrapped in a scikit-learn Pipeline (Imputer → Scaler → SMOTE → Model) and evaluated with **5-Fold Stratified Cross-Validation**. Reported metrics for each fold: Accuracy, F1, Precision, Recall, ROC-AUC. Best model selected by AUC for the fusion step.

N o.	Model	Type	Key Strength	Extra Output
1	Logistic Regression	Linear	Interpretable baseline — coefficient per feature	Feature coefficients
2	Random Forest	Ensemble trees	Robust to outliers — feature importance ranking	Feature importance plot
3	XGBoost	Gradient boosting	Best accuracy on tabular medical data — handles imbalance	SHAP values
4	LightGBM	Gradient boosting	Faster than XGBoost — good for larger filtered datasets	Gain importance
5	SVM (RBF kernel)	Kernel method	Strong on small datasets — effective with scaled features	Support vectors
6	KNN	Distance-based	Non-parametric — useful as reference baseline	Decision boundary
7	MLP Neural Net	Neural network	3 hidden layers (128 → 64 → 32) — learns non-linear patterns	Loss curves

5.2 Cross-Validation Framework

Strategy: 5-Fold Stratified Cross-Validation — each fold preserves the class ratio.

Pipeline order: SimpleImputer → StandardScaler → SMOTE → Model (SMOTE only on train fold)

Metrics per fold: Accuracy, F1 (macro), Precision (macro), Recall (macro), ROC-AUC

Final report: Mean ± Std across 5 folds + test set metrics for the best model

Cross-Validation Code — All 7 Models:

```
import numpy as np from sklearn.linear_model import LogisticRegression from sklearn.ensemble import
RandomForestClassifier from sklearn.svm import SVC from sklearn.neighbors import KNeighborsClassifier from
sklearn.neural_network import MLPClassifier from xgboost import XGBClassifier from lightgbm import
LGBMClassifier from sklearn.impute import SimpleImputer from sklearn.preprocessing import StandardScaler
from sklearn.model_selection import StratifiedKFold, cross_validate from imblearn.over_sampling import
SMOTE from imblearn.pipeline import Pipeline as ImbPipeline # Define all 7 models MODELS = { 'Logistic
Regression': LogisticRegression(max_iter=1000, random_state=42), 'Random Forest':
RandomForestClassifier(n_estimators=200, random_state=42), 'XGBoost': XGBClassifier(n_estimators=200,
eval_metric='logloss', random_state=42), 'LightGBM': LGBMClassifier(n_estimators=200, random_state=42),
'SVM (RBF)': SVC(kernel='rbf', probability=True, random_state=42), 'KNN':
KNeighborsClassifier(n_neighbors=5), 'MLP': MLPClassifier(hidden_layer_sizes=(128,64,32), max_iter=500,
random_state=42), } SCORING = ['accuracy','f1_macro','precision_macro','recall_macro','roc_auc'] cv =
StratifiedKFold(n_splits=5, shuffle=True, random_state=42) results = {} for name, model in MODELS.items():
pipe = ImbPipeline([ ('imputer', SimpleImputer(strategy='median')), ('scaler', StandardScaler()),
('smote', SMOTE(random_state=42)), ('model', model), ]) scores = cross_validate(pipe, X_train, y_train,
cv=cv, scoring=SCORING, return_train_score=True) results[name] = { 'Accuracy':
f"{scores['test_accuracy'].mean():.3f} ± {scores['test_accuracy'].std():.3f}", 'F1':
f"{scores['test_f1_macro'].mean():.3f} ± {scores['test_f1_macro'].std():.3f}", 'AUC':
f"{scores['test_roc_auc'].mean():.3f} ± {scores['test_roc_auc'].std():.3f}", } print(f'{name}: AUC =
{scores["test_roc_auc"].mean():.3f}') # Select best model by AUC best_name = max(results, key=lambda k:
float(results[k]['AUC'].split(' ')[0])) print(f'\nBest model: {best_name}')
```

5.3 Per-Model Evaluation — Additional Outputs

Model	Confusion Matrix	ROC Curve	Learning Curve	Feature Importance
Logistic Regression	Yes	Yes	Yes	Coefficients
Random Forest	Yes	Yes	Yes	Gini importance bar chart
XGBoost	Yes	Yes	Yes	SHAP summary plot
LightGBM	Yes	Yes	Yes	Gain importance
SVM (RBF)	Yes	Yes	Yes	Support vectors only
KNN	Yes	Yes	Yes	K-value sensitivity plot
MLP Neural Net	Yes	Yes	Yes	Train/val loss curves

SECTION 6 — BRANCH B: 8 CNN ARCHITECTURES + TRANSFER LEARNING

6.1 Architecture Comparison

Model	Params	Top-1	Speed	Key Innovation	Best Use
MobileNetV2	3.4 M	72.0%	Fastest	Inverted residuals, depthwise conv	CPU/mobile deployment
AlexNet*	60 M	—	Slow	First deep CNN — trained from scratch	Historical baseline
ResNet50V2	25 M	75.6%	Medium	Skip connections: output = F(x)+x	Robust baseline
InceptionV3	23 M	77.9%	Medium	Parallel 1x1, 3x3, 5x5 convolutions	Multi-scale CXR features
DenseNet121	8 M	74.7%	Fast	All layers connect to all prior layers	CheXNet backbone — recommended
EfficientNetB3	12 M	81.6%	Medium	Compound depth/width/resolution scaling	Best accuracy/size ratio
ViT	86 M	81.8%	Slow	Image patches as transformer tokens	Global bilateral patterns
VGG16	138 M	71.5%	Slowest	Deep uniform 3x3 conv stacks	Heavy reference baseline

* AlexNet has no ImageNet pre-training — trained from scratch on X-ray data only.

6.2 Transfer Learning Strategy — 2-Phase Training

Parameter	Phase 1 — Head Training	Phase 2 — Fine-Tuning
Backbone	ALL layers frozen	Last 20 layers trainable
Learning Rate	1e-3 (higher — head not pre-trained)	1e-5 (very small — prevent forgetting)
Epochs	25 max (early stopping patience=8)	15 max (early stopping patience=6)
Batch Size	32 (GPU) / 16 (CPU)	32 (GPU) / 8 (CPU)
Risk	LR too small → slow convergence	LR too large → catastrophic forgetting

Classification Head — same for all 8 architectures:

```
from tensorflow.keras import Model
from tensorflow.keras.layers import GlobalAveragePooling2D, BatchNormalization, Dense, Dropout, Activation
from tensorflow.keras.applications import DenseNet121 # example
def build_model(base_app, input_shape=(224, 224, 3)):
    base = base_app(include_top=False, weights='imagenet', input_shape=input_shape)
    base.trainable = False # Phase 1: freeze all
    x = GlobalAveragePooling2D()(base.output)
    x = BatchNormalization()(x)
    x = Dense(512, activation='relu')(x)
    x = Dropout(0.5)(x)
    x = Dense(256, activation='relu')(x)
    x = Dropout(0.3)(x) # MUST be float32 - float16 causes NaN loss
    x = Dense(1, dtype='float32')(x)
    output = Activation('sigmoid', dtype='float32')(x)
    return Model(inputs=base.input, outputs=output)
```

SECTION 7 — FUSION STRATEGIES

7.1 The Core Challenge — Datasets Are Not Paired

Key constraint: The two datasets do NOT share the same patients. There is no row in the CSV that corresponds to a specific X-ray image. This means end-to-end joint training is not possible with these datasets. The correct approach is Late Fusion (train separately, combine outputs). Intermediate Fusion requires a paired dataset like MIMIC-CXR.

7.2 Strategy 1 — Late Fusion (Recommended for this project)

Each model is trained completely independently on its own dataset. At inference time, both models produce a pneumonia probability and the outputs are combined with a weighted average. Weights are set to the validation AUC of each model.

```
# Both models are already trained and saved
import numpy as np
def late_fusion_predict(patient_clinical, patient_xray, tabular_model, cnn_model, w1, w2):
    # w1, w2 = validation AUC of each model
    # Branch A - tabular prediction
    p1 = tabular_model.predict_proba(patient_clinical)[0][1]
    # Branch B - CNN prediction
    p2 = float(cnn_model.predict(patient_xray)[0][0])
    # Weighted average fusion
    p_final = (w1 * p1 + w2 * p2) / (w1 + w2)
    # Decision label
    label = 'PNEUMONIA' if p_final > 0.5 else 'NORMAL'
    print(f'Clinical: {p1:.3f} | X-Ray: {p2:.3f} | Fused: {p_final:.3f} → {label}')
    return p_final, label
# Example usage after training:
# w1 = 0.87 (tabular model validation AUC)
# w2 = 0.97 (CNN validation AUC)
# result = late_fusion_predict(X_new, img_new, tab_model, cnn_model, w1, w2)
```

7.3 Strategy 2 — Intermediate Fusion (Advanced — needs paired data)

Instead of combining probabilities, extract the internal embeddings from each model, concatenate them into a joint vector, and train a new classification head on top. The fusion layer learns which cross-modal combinations are most predictive.

```
from tensorflow.keras import Input, Model
from tensorflow.keras.layers import Dense, Dropout, Concatenate, GlobalAveragePooling2D
from tensorflow.keras.applications import EfficientNetB3

# BRANCH A - Tabular MLP producing 64-dim embedding
tabular_input = Input(shape=(16,), name='clinical_input')
x1 = Dense(128, activation='relu')(tabular_input)
x1 = Dropout(0.3)(x1)
tabular_emb = Dense(64, name='clinical_embedding')(x1)

# [64] # BRANCH B - CNN producing 64-dim embedding
base_cnn = EfficientNetB3(include_top=False, weights='imagenet', input_shape=(224, 224, 3))
x2 = GlobalAveragePooling2D()(base_cnn.output)
x2 = Dense(128, activation='relu')(x2)
image_emb = Dense(64, name='image_embedding')(x2)

# [64] # FUSION - concatenate: [64] + [64] = [128]
merged = Concatenate()([tabular_emb, image_emb])
merged = Dense(64, activation='relu')(merged)
merged = Dropout(0.4)(merged)
output = Dense(1, activation='sigmoid', dtype='float32')(merged)

# Joint model - needs PAIRED data (same patient in both modalities)
fusion_model = Model(inputs=[tabular_input, base_cnn.input], outputs=output, name='intermediate_fusion')

# For unpaired data: use Late Fusion above.
# For paired data (MIMIC-CXR): use this intermediate fusion.
```

Property	Late Fusion	Intermediate Fusion
What is combined	Output probabilities p1 and p2	Internal embedding vectors (64-dim each)
Data required	Independent datasets — no pairing needed	Paired: same patient in both CSV and X-ray
Learns cross-modal	No — static weighted average	Yes — fusion layer learns interactions
Complexity	Low — 1 line of code	High — joint training required
For this project	RECOMMENDED — use now	Future work with MIMIC-CXR

SECTION 8 — EVALUATION METRICS & VALIDATION

8.1 Why Accuracy Alone Is Not Enough

In pneumonia diagnosis, a missed case (false negative) is far more dangerous than a false alarm (false positive). Standard accuracy is misleading with class imbalance and does not distinguish these clinically different errors. A full set of metrics is required.

Metric	Formula	Clinical Meaning	Target
Accuracy	$(TP+TN)/(TP+TN+FP+FN)$	Overall correct classification rate	>90%
Sensitivity/Recall	$TP/(TP+FN)$	Rate of real pneumonia caught — MOST CRITICAL	>95%
Specificity	$TN/(TN+FP)$	Rate of normal cases correctly identified	>85%
Precision	$TP/(TP+FP)$	% of predicted pneumonia that are truly pneumonia	>85%
F1 Score	$2 \times (P \times R) / (P + R)$	Harmonic balance of precision and recall	>90%
ROC-AUC	Area under ROC curve	Overall discriminative power — primary metric	>0.95
PR-AUC	Area under PR curve	More reliable than ROC under class imbalance	>0.90

8.2 12 Plots Generated Per Model

No	Plot	What It Shows
1	Histograms + KDE	Distribution of each clinical feature — skew, outliers, normality
2	Boxplots by label	Feature values split by Normal vs Pneumonia — t-test p-value shown
3	Violin plots	Full distribution shape and density by diagnosis label
4	Target balance pie + bar	Label balance, severity breakdown, sex distribution
5	Binary signs bar	Prevalence of each binary clinical sign (cough, fever, etc.)
6	Correlation heatmap	Full Pearson correlation matrix — feature interdependencies
7	Missing values heatmap	Pattern and percentage of missing values per feature
8	Age analysis KDE + box	Age distribution colored by severity level
9	Pairplot (key vitals)	Scatter matrix of spo2, temperature, heart_rate colored by label
10	Confusion matrix	TP/TN/FP/FN per model on test set
11	ROC + PR curves (all 7)	Combined ROC and Precision-Recall curves for all models on one plot
12	Feature importance / SHAP	Which clinical features most influence each model's prediction

SECTION 9 — IMPLEMENTATION ROADMAP & CURRENT STATUS

Phase	Name	Status	Key Deliverables
1	Literature Review	COMPLETE	CNN architectures, tabular ML, fusion strategies, medical AI papers studied
2	EDA — Both Datasets	COMPLETE	Age distributions, class imbalance identified, validation set problem found
3	CSV Pediatric Filter	COMPLETE	Filter to ages 1–5 coded and verified. pediatric_pneumonia.csv created.
4	Preprocessing — Both Branches	COMPLETE	Branch A: encode+scale+SMOTE pipeline. Branch B: resize+normalize+augment+fix split.
5	Branch A Training (7 models)	NOT STARTED	Run 5-fold CV for all 7 tabular models. Select best by AUC.
6	Branch B Training (8 CNNs)	NOT STARTED	Train 8 CNN architectures with 2-phase transfer learning. Select best by val AUC.
7	Late Fusion Implementation	DESIGNED	Code ready — awaits trained models. Weight by AUC of each branch.
8	Full Evaluation + Report	PLANNED	Confusion matrices, ROC curves, SHAP, final comparison of all models.
9	Multi-Dataset Expansion	FUTURE	Add NIH ChestX-ray14 + CheXpert for generalised model beyond pediatric.
10	Intermediate Fusion (paired data)	FUTURE	MIMIC-CXR paired data. Joint embedding training.

SECTION 10 — CONCLUSION & RECOMMENDATIONS

This report has described the complete design of a multimodal AI pneumonia detection system combining clinical tabular data and chest X-ray images. The dual-branch architecture mirrors real clinical practice and addresses a critical diagnostic need in under-resourced settings like Algeria, where specialist radiologists are scarce and diagnostic delays are long.

Topic	Conclusion
CSV age problem	Filter to ages 1–5 before training. Use pneumonia_dataset_pediatric.csv for Branch A only.
Best fusion now	Late Fusion is the correct strategy for these unpaired datasets. Simple and effective.
Best tabular model (expected)	XGBoost or LightGBM based on clinical tabular literature. Confirmed after 5-fold CV.
Best CNN (expected)	DenseNet121 (CheXNet backbone) or EfficientNetB3. Confirmed after training.
Current status	Literature review and preprocessing complete. Branch A and B training is the next step.
Future novelty	Confidence-aware rejection (uncertain cases sent to radiologist), severity scoring (0–3), multi-dataset expansion.
Publication potential	The Algerian-context multimodal system with local EMR integration is a genuine contribution. Publishable in IEEE Access or Frontiers in Digital Health after training results.

Contact for more information, collaboration, or questions:

Kassouar Fatima (Project Manager) — kassouarfatima48@gmail.com
GitHub: github.com/AminaMar/pediatric-pneumonia-detection
University of Saida, Algeria · 2025–2026